



Ambulance Stretcher come with fixed devices

PM-3D-WC

It made of high-strength aluminium alloy w/Waterproof foamed cushion.
The mechanical folding legs are controlled by handles on bothside.
Can be arbitrarily adjusted by gas spring and it can be converted into a wheel chair.

Key Specifications :

Unfolded Size : 191(L) x 53(W) x 82(H) cm

Adjustable Level : 2 high level
Wheel chair mode

Folded Size : 191(L) x 53(W) x 25(H) cm

Self Weight : 39 kg

Load Bearing : 159 kg



2 adjustable high level



Wheel chair mode

