



Ambulance Stretcher come with fixed devices

PM-3H-WF



It made of high-strength aluminium alloy w/Waterproof foamed cushion.
Can be adjusted six different high position & surface can be used as carrying stretcher when it separated from main part.
The mechanical folding legs are controlled by handles on both side. Only one person can operate the stretcher when sending it to the ambulance.

Key Specifications :

Unfolded Size : 191(L) x 53(W) x 82(H) cm

Adjustable Level : 4 high level
Removable mode
Downhill mode

Folded Size: : 191(L) x 53(W) x 25(H) cm

Self Weight : 39 kg

Load Bearing : 159 kg

4 adjustable high level

