



First Aid Kit Metal Medium Casing

PM-01-MM



Key Specifications :

Size : 290(L) x 90(W) x 135(H) mm

Weight : 1.6 kg (with Content)

Total Listing Item : 32 list item

Content Listing :

1.	Quality Steel Container		1 pc		
	CONTENT				
2.	Sterile Finger Dressing No.7	5cm x 5cm	2 pcs		
3.	Sterile Lint Dressing No. 8	8cm x 9cm	1 pc		
4.	Sterile Eye Pad	6cm x 8cm	1 pc		
5.	Confirming Elastic Gauze Bandage	6cm x 4m	1 roll		
6.	W.O.W. Bandage	2.5cm x 5m	1 roll		
7.	W.O.W. Bandage	5cm x 5m	1 roll		
8.	W.O.W. Bandage	7.5cm x 5m	1 roll		
9.	Cotton Balls 10's		1 box		
10.	Gauze Swabs 10's	5cm x 5cm	1 box		
11.	Cotton Buds 100's		1 pkt		
12.	Triangular Bandage	91.5cm x 91.5cm x 130cm	1 pc		
13.	Sterile Assorted Plaster Strip 20's/box		1 set		
	8 x First Aid Bandage 19mm x 72mm				
	4 x First Aid Bandage 25mm x 72mm				
	4 x First Aid Bandage 38mm x 38mm (Square)				
	2 x First Aid Bandage 45mm x 55mm (Butterfly Shape)				
	2 x First Aid Bandage 38mm x 76mm (H Shape)				
	Assorted Cleaning Wipe				
14.	Alcohol Wipe	30mm x 65mm	4 pcs		
15.	Normal Saline Wipe	30mm x 65mm	2 pcs		
16.	Sting Relief Wipe	30mm x 65mm	1 pc		
17.	Antiseptic Cream	15 gm	1 pc		
18.	Adhesive Tape w/dispenser	1.25cm x 5m	1 pc		
19.	S/Steel Utility Scissors Curve	15cm	1 pc		
20.	Plastic Tweezers	12cm	1 pc		
21.	Assorted Safety Pin 12's				1 set
22.	Bio Hazard Zip Lock Bag	8" x 11"			1 pc
23.	Vinyl Gloves 2's				1 pair
24.	General Usage Card				1 pc
25.	Basic Life Support Card				1 pc
26.	CPR & First Aid Guide Booklet				1 pc
27.	Nail				2 pcs
28.	Content Listing Leaflet				1 pc