



Four Fold Stretcher (w/sewn 2 strap)

PM-1A9-FF



Frame made of high-strength aluminium alloy w/Durable Waterproof PVC Coated Nylon Canvas 1000D Orange color w/sewn 2 strap in carrying bag. It's light-weighted, compact and totally portable.

Key Specifications :

Unfolded Size :	221(L) x 53(W) x 15(H) cm
Folded Size :	61(L) x 30(W) x 12(H) cm
Self Weight :	6.2 kg
Load Bearing :	159 kg