



Hospital Mechanical Hilo Bed (Single Fowler)



BA-3010-PN

- Single Fowler *C/W PE Blow head & foot board
- Steel Construction
- With adjustable Height (HI-LO Position) Backrest & Kneerest by mean of 3 Crank Winding mechanism at foot end of bed
- Mobile on 4 x 125mm wheel Diagonally fitted w/Brakes
- Provision for IV Drip pole attachment
- Finishing of hard epoxy powder coating
- HILO Position from 420mm - 720mm (w/o mattress)
- Mattress Size : L 1905mm x W 864 mm * DIM : LI 2090 x W 955 x H 420-740 mm

Key Specifications :

Hilo position from floor to mattress support top : 420-740mm

Backrest : 0° - 80°

Size : 211(L) x 95.5(W) x 57(H) cm

